

CD167a(DDR1)

产品信息

产品名称	产品货号	产品规格
CD167a(DDR1)	ORS7227PAb	50 μ l/100 μ l/1ml

产品简介

CD167a is a membrane type II receptor kinase, containing a factor VIII-like domain. At least five isoforms are generated through alternative splicing of the kinase domain of human DDR1 gene. These isoforms are named with suffixes a to e. CD167a is expressed on epithelial cells, keratinocytes, leukocytes, monocytes and has been reported to be overexpressed in some breast carcinomas. CD167a expression can be upregulated by p53. CD167a has been reported to interact with a variety of proteins including collagen type II alpha 1, collagen type III alpha 1, collagen type V alpha 2, collagen type XI alpha 1, phospholipase gamma 1, SHC, and the lipid-anchored docking protein FRS2. CD167a can be phosphorylated on Y513. Animal studies suggest that CD167a is involved in arterial wound healing.

产品性质

Alias: Discoidin domain receptor, DDR1, neuroepithelial tyrosine kinase, epithelial-specific receptor kinase, tyrosine kinase receptor E (trkE), cell adhesion kinase (CAK), PTK3, RTK6

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q08345

Clone: 6E5

Entrez GeneID: 780

Host\Isotype: Mouse IgG2a

Reactivity: Human

Purification: affinity purified by Protein G

Applications: FCM=1/100-1/200

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信