

## CD20

### 产品信息

| 产品名称 | 产品货号       | 产品规格                       |
|------|------------|----------------------------|
| CD20 | ORS7230PAb | 50 $\mu$ l/100 $\mu$ l/1ml |

### 产品简介

CD20 is a 33-37 kD protein, a member of the MS4A family, with four transmembrane spanning regions that present as a homo-oligomeric complexes in the cell surface when associating with MHC class I and II, CD53, CD81, and CD82. CD20 is expressed on B cells and a subset of T cells, but not on plasma cells. CD20 regulates B-cell activation and proliferation. Its ligation promotes transmembrane  $Ca^{2+}$  trafficking. CD20 is an important therapeutic target in the treatment of B cell lymphomas and leukemias.

### 产品性质

**Alias:** B1, Bp35, MS4A1, CVID5, FMC7, LEU-16, S7

**Clonality:** Recombinant Mouse Monoclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** P11836

**Clone:** 4H1

**Entrez GeneID:** 931

**Host\Isotype:** Mouse IgG2a

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** FCM=1/100-1/200

**Immunogen:** -

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信