

FLAG-tag Antibody-AF633

产品信息

产品名称	产品货号	产品规格
FLAG-tag Antibody-AF633	ORS7182PAb	100 μ l/500 μ l/1ml

产品简介

Mouse monoclonal antibody to FLAG-tag-AF633 labeled, The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads., Recognizes C-terminal, internal, and N-terminal FLAG-tag fusion proteins.

产品性质

Alias: Anti-FLAG-tag Antibody-AF633 labeled; Anti FLAG tag Antibody AF633 labeled

Clonality: Monoclonal Antibody

Antibody Type: Secondary Antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: Mouse

Reactivity: -

Purification: affinity purified by Protein G

Applications: WB=1/2000 - 1/5000; ICC=1/500 - 1/2000; FCM=1/1000 - 1/4000

Immunogen: KLH-conjugated synthetic peptide encompassing a sequence of FLAG-tag.

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信