

Flag Tag

产品信息

产品名称	产品货号	产品规格
Flag Tag	ORS0120PAb	100 μ l/50 μ l

产品简介

The antibody is useful for the identification of proteins with Flag markers at free N-termini and C-termini by common immunological applications in yeast, animal, and E. coli cells.

产品性质

Alias: -

Clonality: Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: 2F11G1

Entrez GeneID: -

Host\Isotype: Mouse IgG2b

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/500-1/1000; IP=1/10-1/20; IHC=1/100 - 1/1000; ICC=1/100-1/200; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Synthetic peptide corresponding to aa(DYKDDDDK), conjugated to KLH.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信