

HSP60

产品信息

产品名称	产品货号	产品规格
HSP60	ORS0417PAb	100 μ l/50 μ l

产品简介

This gene encodes a member of the chaperonin family. The encoded mitochondrial protein may function as a signaling molecule in the innate immune system. This protein is essential for the folding and assembly of newly imported proteins in the mitochondria. This gene is adjacent to a related family member and the region between the 2 genes functions as a bidirectional promoter. Several pseudogenes have been associated with this gene. Two transcript variants encoding the same protein have been identified for this gene. Mutations associated with this gene cause autosomal recessive spastic paraplegia 13. (provided by RefSeq)

产品性质

Alias: HLD4; CPN60; GROEL; HSP60; HSP65; SPG13; HSP-60; HuCHA60; HSPD1

Clonality: Monoclonal Antibody

Antibody Type: Loading Control Antibody

Swissprot: P10809

Clone: 3G8

Entrez GeneID: 3329

Host\Isotype: Mouse IgG1

Reactivity: Human, Mouse, Rat, Monkey

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IHC=1/200 - 1/1000; ICC=1/200 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fragment of human HSP60 expressed in E. Coli.

Formulation: Ascitic fluid containing 0.03% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信