

Human IgA - PE

产品信息

产品名称	产品货号	产品规格
Human IgA - PE	ORS7115PAb	100 μ l/500 μ l/1ml

产品简介

Mouse Monoclonal Secondary Antibody to Human IgA PE labeled, The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads., By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgA. No antibody was detected against non immunoglobulin serum proteins.

产品性质

Alias: Mouse Anti-Human IgA - PE; Mouse Anti Human IgA PE

Clonality: Monoclonal Antibody

Antibody Type: Secondary Antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: Mouse

Reactivity: Human

Purification: affinity purified by Protein G

Applications: IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/1000 - 1/4000; ELISA=Use at an assay dependent concentration

Immunogen: Human IgA

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信