

## Human IgG Fab - APC

### 产品信息

产品名称	产品货号	产品规格
Human IgG Fab - APC	ORS7122PAb	100 $\mu$ l/500 $\mu$ l/1ml

### 产品简介

Mouse Monoclonal Secondary Antibody to Human IgG Fab APC labeled, The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads., By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fab. No antibody was detected against non immunoglobulin serum proteins.

### 产品性质

**Alias:** Mouse Anti-Human IgG Fab - APC; Mouse Anti Human IgG Fab APC

**Clonality:** Monoclonal Antibody

**Antibody Type:** Secondary Antibody

**Swissprot:** -

**Clone:** -

**Entrez GeneID:** -

**Host\Isotype:** Mouse

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/1000 - 1/4000; ELISA=Use at an assay dependent concentration

**Immunogen:** Human IgG Fab

**Formulation:** Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信