

MBP Tag

产品信息

产品名称	产品货号	产品规格
MBP Tag	ORS0128PAb	100 μ l/50 μ l

产品简介

Epitope tagging offers an easy and universal strategy for the identification and purification of proteins derived by recombinant DNA technology. The insertion of a Maltose Binding Protein (MBP) tag creates a stable fusion product that does not interfere with the bioactivity of the protein or with the biodistribution of the MBP tagged product. Cleavage by factor Xa separates MBP from its partner protein.

产品性质

Alias: Mal

Clonality: Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: 3G1A3

Entrez GeneID: -

Host\Isotype: Mouse IgG1

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/500-1/1000; IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Recombinant fusion protein with Maltose binding protein tag.

Formulation: Ascitic fluid containing 0.03% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信