

Metapneumovirus

产品信息

产品名称	产品货号	产品规格
Metapneumovirus	ORS4834PAb	1mg

产品简介

-

产品性质

Alias: -

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: -

Reactivity: Matrix Protein (genotypes A and B)

Purification: affinity purified by Protein G

Applications: -

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信