

## PCNA

### 产品信息

产品名称	产品货号	产品规格
PCNA	ORS2443PAb	100 $\mu$ l/50 $\mu$ l

### 产品简介

The protein encoded by this gene is found in the nucleus and is a cofactor of DNA polymerase delta. The encoded protein acts as a homotrimer and helps increase the processivity of leading strand synthesis during DNA replication. In response to DNA damage, this protein is ubiquitinated and is involved in the RAD6-dependent DNA repair pathway. Two transcript variants encoding the same protein have been found for this gene. Pseudogenes of this gene have been described on chromosome 4 and on the X chromosome. <br /><br />

### 产品性质

Alias: -

Clonality: Monoclonal Antibody

Antibody Type: Loading Control Antibody

Swissprot: P12004

Clone: 7H4F8

Entrez GeneID: 5111

Host\Isotype: Mouse IgG1

Reactivity: Human, Mouse, Rat, Rabbit, Monkey

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IHC=1/100 - 1/500; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fragment of human PCNA (AA: 53-196 ) expressed in E. Coli.

Formulation: Ascitic fluid containing 0.03% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>  
<https://www.aoruicell.com>



手机号: 13061929181(微信)  
 13061929770(微信)



邮箱: sale4@aoruisai.com  
 sale5@aoruisai.com



电话: 021-57896695  
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信