

TUBB1

产品信息

产品名称	产品货号	产品规格
TUBB1	ORS2141PAb	100 μl/50 μl

产品简介

This gene encodes a member of the beta tubulin protein family. Beta tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. This protein is specifically expressed in platelets and megakaryocytes and may be involved in proplatelet production and platelet release. A mutations in this gene is associated with autosomal dominant macrothrombocytopenia. Two pseudogenes of this gene are found on chromosome Y.

产品性质

Alias: M40; TUBB1; TUBB5; CDCBM6; CSCSC1; OK/SW-cl.56

Clonality: Monoclonal Antibody

Antibody Type: Loading Control Antibody

Swissprot: P07437

Clone: 4D4G6

Entrez GeneID: 203068

Host\Isotype: Mouse IgG1

Reactivity: Human, Mouse, Monkey

Purification: affinity purified by Protein G

Applications: WB=1/1000 - 1/2000; IHC=1/200 - 1/1000; ICC=1/200 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fragment of human TUBB1 (AA: 1-444) expressed in E. Coli.

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信