

beta-Actin

产品信息

产品名称	产品货号	产品规格
beta-Actin	ORS0146PAb	100 μ l/50 μ l

产品简介

Beta-actin (PS1TP5-binding protein 1), also known as ACTB, PS1TP5BP1. Entrez Protein NP_001092. It is one of six different actin proteins. Actin, a ubiquitous eukaryotic protein, is the major component of the cytoskeleton. Actins are highly conserved proteins that are involved in various types of cell motility, structure, and integrity. Actin is ubiquitously expressed in all eukaryotic cells. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

产品性质

Alias: PS1TP5BP1; ACTB

Clonality: Monoclonal Antibody

Antibody Type: Loading Control Antibody

Swissprot: P60709

Clone: 8H10D10

Entrez GeneID: 60

Host\Isotype: Mouse IgG2b

Reactivity: Human, Mouse, Hamster, Rat, Monkey, Dog

Purification: affinity purified by Protein G

Applications: WB=1/10000 - 1/40000; IP=1/100-1/200; ICC=1/50 - 1/500; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Synthetic peptide corresponding to amino-terminal residues of human beta-Actin, conjugated to KLH.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信