

ADH6

产品信息

产品名称	产品货号	产品规格
ADH6	ORS8633MAb	100 μ l/50 μ l

产品简介

ADH6 encodes class V alcohol dehydrogenase, which is a member of the alcohol dehydrogenase family. Members of this family metabolize a wide variety of substrates, including ethanol, retinol, other aliphatic alcohols, hydroxysteroids, and lipid peroxidation products. This protein is expressed in the stomach as well as in the liver, and it contains a glucocorticoid response element upstream of its 5' UTR, which is a steroid hormone receptor binding site.

产品性质

Alias: Alcohol dehydrogenase 6, ADH6

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P28332

Clone: -

Entrez GeneID: 130

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000; IHC=1/100-1/500; FCM=1/10-1/50

Immunogen: This ADH6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 209-240 amino acids from the Central region of human ADH6.

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol, prepared by Saturated Ammonium Sulfate (SAS) .

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信