

CDC73

产品信息

产品名称	产品货号	产品规格
CDC73	ORS8628MAb	100 μ l/50 μ l

产品简介

-

产品性质

Alias: Parafibromin, Cell division cycle protein 73 homolog, Hyperparathyroidism 2 protein, CDC73, Clorf28, HRPT2

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q6P1J9

Clone: -

Entrez GeneID: -

Host\Isotype: -

Reactivity: Human, Mouse, Rat

Purification: affinity purified by Protein G

Applications: -

Immunogen: This CDC73 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 132-161 amino acids from the Central region of human CDC73.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信