

CEP63

产品信息

产品名称	产品货号	产品规格
CEP63	ORS8720MAb	100 μ l/50 μ l

产品简介

CEP63 is a protein with six coiled-coil domains. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells.

产品性质

Alias: Centrosomal protein of 63 kDa, Cep63, CEP63

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q96MT8

Clone: -

Entrez GeneID: 80254

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000; IHC=1/100-1/500

Immunogen: This CEP63 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 520-548 amino acids from the Central region of human CEP63.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信