

澳睿赛---专业多元化试剂厂商

AoRuiSai----Professional and diversified reagent manufacturers

ELISA试剂盒/生化检测试剂盒/重组蛋白/抗体等

CLM1 (N-term)

产品信息

产品名称	产品货号	产品规格
CLM1 (N-term)	ORS8742MAb	100 μ 1/50 μ 1

产品简介

CLM1 is an inhibitory receptor of the Ig superfamily expressed on myeloid cells. It mediates negative regulatory signals by recruiting SHP1 (PTPN6; MIM 176883) or SHIP.

产品性质

Alias: CMRF35-like molecule 1, CLM-1, CD300 antigen-like family member F, Immune receptor expressed on myeloid cells 1, IREM-1, Immunoglobulin superfamily member 13, IgSF13, NK inhibitory receptor, CD300f, CD300LF, CD300F, CLM1, IGSF13, IREM1, NKIR

Clonality: Polyclonal Antibody Antibody Type: Primary antibody

Swissprot: Q8TDQ1

Clone: -

Entrez GeneID: 146722 Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G Applications: WB=1/1000; IHC=1/100-1/500

Immunogen: This CLM1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 56-84 amino acids from the N-terminal region of human CLM1

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid

freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



澳睿赛---专业多元化试剂厂商

AoRuiSai----Professional and diversified reagent manufacturers

ELISA试剂盒/生化检测试剂盒/重组蛋白/抗体等

- 2. Avoid repeated freeze-thaw cycles.
- 3. For your safety and health, please wear a lab coat and disposable gloves



● 网站: https://www.aoruisai.com ☐手机号: 13061929181(微信)

https://www.aoruicell.com 13061929770(微信)

021-57875552 务/监督投诉电话: 13122716991

