

## DDO

### 产品信息

| 产品名称 | 产品货号       | 产品规格                   |
|------|------------|------------------------|
| DDO  | ORS8718MAb | 100 $\mu$ l/50 $\mu$ l |

### 产品简介

DDO is a peroxisomal flavoprotein that catalyzes the oxidative deamination of D-aspartate and N-methyl D-aspartate. Flavin adenine dinucleotide or 6-hydroxyflavin adenine dinucleotide can serve as the cofactor in this reaction.

### 产品性质

**Alias:** D-aspartate oxidase, DASOX, DDO, DDO

**Clonality:** Polyclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** Q99489

**Clone:** -

**Entrez GeneID:** 8528

**Host\Isotype:** Rabbit IgG

**Reactivity:** Human, Mouse

**Purification:** affinity purified by Protein G

**Applications:** WB=1/1000; IHC=1/100-1/500

**Immunogen:** This DDO antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 102-130 amino acids from the Central region of human DDO.

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)



邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)



电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信