

## FBX03

### 产品信息

| 产品名称  | 产品货号       | 产品规格                   |
|-------|------------|------------------------|
| FBX03 | ORS8630MAb | 100 $\mu$ l/50 $\mu$ l |

### 产品简介

FBX03 encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs.

### 产品性质

**Alias:** F-box only protein 3, FBX03, FBX3

**Clonality:** Polyclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** Q9UK99

**Clone:** -

**Entrez GeneID:** 26273

**Host\Isotype:** Rabbit IgG

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** WB=1/1000;FCM=1/10-1/50

**Immunogen:** This FBX03 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 395-422 amino acids from the C-terminal region of human FBX03.

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信