

FLAD1

产品信息

产品名称	产品货号	产品规格
FLAD1	ORS8631MAb	100 μ l/50 μ l

产品简介

FLAD1 encodes the enzyme that catalyzes adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme.

产品性质

Alias: FAD synthase, FAD pyrophosphorylase, FMN adenylyltransferase, Flavin adenine dinucleotide synthase, Molybdenum cofactor biosynthesis protein-like region, FAD synthase region, FLAD1

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q8NFF5

Clone: -

Entrez GeneID: 80308

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000

Immunogen: This FLAD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 551-580 amino acids from the C-terminal region of human FLAD1.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信