

KIR3DL3

产品信息

产品名称	产品货号	产品规格
KIR3DL3	ORS8625MAb	100 μ l/50 μ l

产品简介

KIR3DL3 is transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. This protein are polymorphic and highly homologous and they are found in a cluster on chromosome 19q13.4 within the 1 Mb leukocyte receptor complex (LRC).

产品性质

Alias: Killer cell immunoglobulin-like receptor 3DL3, CD158 antigen-like family member Z, Killer cell inhibitory receptor 1, CD158z, KIR3DL3, CD158Z, KIR3DL7, KIRC1

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q8N743

Clone: -

Entrez GeneID: -

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000;IHC=1/100-1/500;FCM=1/10-1/50

Immunogen: This KIR3DL3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 132-158 amino acids from the Central region of human KIR3DL3.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信