

NARFL

产品信息

产品名称	产品货号	产品规格
NARFL	ORS8706MAb	100 μ l/50 μ l

产品简介

—

产品性质

Alias: Cytosolic Fe-S cluster assembly factor NARFL, Iron-only hydrogenase-like protein 1, IOP1, Nuclear prelamins A recognition factor-like protein, Protein related to Narf, NARFL, PRN

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q9H6Q4

Clone: —

Entrez GeneID: 64428

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000;FCM=1/10-1/50

Immunogen: This NARFL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 240-269 amino acids from the Central region of human NARFL.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: —

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信