

PXMP3

产品信息

产品名称	产品货号	产品规格
PXMP3	ORS8627MAb	100 μ l/50 μ l

产品简介

PXMP3 encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import.

产品性质

Alias: Peroxisome biogenesis factor 2, 35 kDa peroxisomal membrane protein, Peroxin-2, Peroxisomal membrane protein 3, Peroxisome assembly factor 1, PAF-1, RING finger protein 72, PEX2, PAF1, PMP3, PMP35, PXMP3, RNF72

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P28328

Clone: -

Entrez GeneID: 5828

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse

Purification: affinity purified by Protein G

Applications: WB=1/1000; IHC=1/100-1/500; FCM=1/10-1/50

Immunogen: This PXMP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 165-194 amino acids from the Central region of human PXMP3.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信