

RPS13

产品信息

产品名称	产品货号	产品规格
RPS13	ORS8719MAb	100 μ l/50 μ l

产品简介

RPS13, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S15P family of ribosomal proteins. It is located in the cytoplasm. The protein has been shown to bind to the 5.8S rRNA in rat. The gene product of the E. coli ortholog (ribosomal protein S15) functions at early steps in ribosome assembly. This gene is co-transcribed with two U14 small nucleolar RNA genes, which are located in its third and fifth introns. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

产品性质

Alias: 40S ribosomal protein S13, RPS13

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P62277

Clone: -

Entrez GeneID: 6207

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000; IHC=1/100-1/500; FCM=1/10-1/50

Immunogen: This RPS13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 82-110 amino acids from the Central region of human RPS13.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.



网站: <https://www.aoruisai.com>

<https://www.aoruicell.com>



手机号: 13061929181(微信)

13061929770(微信)



邮箱: sale4@aoruisai.com

sale5@aoruisai.com



电话: 021-57896695

021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信