

SPRN (C-term)

产品信息

产品名称	产品货号	产品规格
SPRN (C-term)	ORS8713MAb	100 μ l/50 μ l

产品简介

SPRN has PrP(C)-like neuroprotective activity. SPRN may act as a modulator for the biological actions of normal and abnormal PrP.

产品性质

Alias: Shadow of prion protein, Protein shadoo, SPRN, SHO

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q5BIV9

Clone: -

Entrez GeneID: 503542

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/1000

Immunogen: This SPRN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 76-105 amino acids from the C-terminal region of human SPRN.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信