

## TBCD7

### 产品信息

产品名称	产品货号	产品规格
TBCD7	ORS8736MAb	100 $\mu$ l/50 $\mu$ l

### 产品简介

TBC1D7 belongs to the family of proteins sharing a 180- to 200-amino acid TBC domain presumed to have a role in regulating cell growth and differentiation. These proteins share significant homology with TRE2.

### 产品性质

**Alias:** TBC1 domain family member 7, Cell migration-inducing protein 23, TBC1D7, TBC7

**Clonality:** Polyclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** Q9PON9

**Clone:** -

**Entrez GeneID:** 107080638;51256

**Host\Isotype:** Rabbit IgG

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** WB=1/1000

**Immunogen:** This TBCD7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 150-179 amino acids from the Central region of human TBCD7.

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)



邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)



电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信