

TOR1B

产品信息

产品名称	产品货号	产品规格
TOR1B	ORS8716MAb	100 μ l/50 μ l

产品简介

TOR1B may serve as a molecular chaperone assisting in the proper folding of secreted and/or membrane proteins.

产品性质

Alias: Torsin-1B, Torsin ATPase-1B, 364-, Torsin family 1 member B, TOR1B, DQ1

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: 014657

Clone: -

Entrez GeneID: 27348

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse

Purification: affinity purified by Protein G

Applications: WB=1/1000

Immunogen: This TOR1B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 210-237 amino acids from the Central region of human TOR1B.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信