

TSPAN6

产品信息

产品名称	产品货号	产品规格
TSPAN6	ORS8636MAb	100 μ l/50 μ l

产品简介

The protein is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein and is highly similar in sequence to the transmembrane 4 superfamily member 2.

产品性质

Alias: Tetraspanin-6, Tspan-6, A15 homolog, Putative NF-kappa-B-activating protein 321, T245 protein, Tetraspanin TM4-D, Transmembrane 4 superfamily member 6, TSPAN6, TM4SF6

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: 043657

Clone: -

Entrez GeneID: 7105

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse

Purification: affinity purified by Protein G

Applications: WB=1/1000;FCM=1/10-1/50

Immunogen: This TSPAN6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 176~205 amino acids from the C-terminal region of human TSPAN6.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.



网站: <https://www.aoruisai.com>

<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信