

ZRANB1

产品信息

产品名称	产品货号	产品规格
ZRANB1	ORS8707MAb	100 μ l/50 μ l

产品简介

ZRANB1 is a positive regulator of the Wnt signaling pathway that specifically cleaves 'Lys-63'-linked ubiquitin chains. It acts by deubiquitinating APC protein, a negative regulator of Wnt-mediated transcription and may also modulate TNF- α signaling.

产品性质

Alias: Ubiquitin thioesterase ZRANB1, TRAF-binding domain-containing protein, hTrabid, Zinc finger Ran-binding domain-containing protein 1, ZRANB1, TRABID

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q9UG10

Clone: -

Entrez GeneID: 54764

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse

Purification: affinity purified by Protein G

Applications: WB=1/1000; IHC=1/100-1/500

Immunogen: This ZRANB1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 665-692 amino acids from the C-terminal region of human ZRANB1.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信