

# 重组蛋白

Catalog Number:ORS5416NSh

Sheep Immunoglobulin G

Organism Species: Purified from pooled normal sheep  
serum



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

## 【Product Information】

Gene Names: -

Uniprot:-

Purity: >95%(SDS-PAGE)

Endotoxin: -

Source: -

Protein Length: -

Expression Region: -

Sequence: Full length protein

Tag Info:

Form:Liquid

Buffer formulation:PBS, pH 7.4

Applications:ELISA, WB, SDS-PAGE

Predicted Molecular Mass:-

Accurate Molecular Mass:-

Predicted isoelectric point:-

## 【Storage and Stability】

**Storage:** Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of


lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

## 【Identification】


 网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991





### 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

