

# 重组蛋白

Catalog Number:ORSBP0287R0

Recombinant Culex quinquefasciatus Submitted name:

Odorant binding protein OBP28(6050687)

Organism Species: Culex quinquefasciatus (Southern house mosquito) (Culex pungens)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



**【Product Information】**

Gene Names: 6050687

Uniprot: B0XC79

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: Baculovirus

Protein Length: Full Length of Mature Protein

Expression Region: 21-150aa

Sequence:

DEASDKEQAKEQAKQMLRSMTQKCKEAE GASDDDVEAMIDDVMPESQVQKCFHSCVQQQF  
GVSDGQKFLQQGFLEIMMAVGNDEQQQGHAKVAEECDGVANEDRCQLAVDIMTCVKQG  
MEKRGMKVDR

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form: Liquid or Lyophilized powder

Buffer formulation: -

Applications: -

Predicted Molecular Mass: -

Accurate Molecular Mass: 18.6 kDa

Predicted isoelectric point: -

**【Storage and Stability】**

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)


邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



## 【Identification】



## 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

 网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

