

# 重组蛋白

Catalog Number:ORSBP2895R0

Recombinant *Saccharomyces cerevisiae* Alpha-1,3-  
mannosyltransferase MNT3 (MNT3), partial

Organism Species: *Saccharomyces cerevisiae* (strain ATCC  
204508 / S288c) (Baker's yeast)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

## 【Product Information】

Gene Names: MNT3

Uniprot:P40549

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: Baculovirus

Protein Length: Partial

Expression Region: 32-630aa

Sequence:

SRNVTSSNKLNNHASERTAVESSAFNWIEKRQHQVRSENLMNRLSAYFLPFLSRSSHKER  
VLLRQLGNNEIAKSDKCRYIFEVLYKIDPDWDNAQTAKFYNVVDGVDNTLASLLGERLSY  
DYCFLSGQLDPTAIFANSTVNPDLQNRMFPLKKEINESKTVMWPIITDMMTGEAVPAP  
EVDMESSNFNGFWSNWNRLSKGRGFVLTIAEKDVPLFLKQLKVMFESKNELEPFQIVSTG  
NELSAESIAKISSETAKETEQRVYLVCSTVLDTNFANTYISFFQNKWVATLNFNTFEEYIL  
LDADVVPFVGSDFYFFDPSYRESGILLFKDRVMENEQTFQYCIEMLENEVEPSAQERRFIG  
SRLVFDSSLPFSSETSEEASVYYNFFKRLRLHHVDSGLVVVNKLEKLNGLLMSFMLNLDG  
KLQRCVYGDKEIFWLGQLYAGQDYSINPVDGSIIGPVNEEPENDDGHKSGMYIICSTQIA  
HSDSKNRLWVNGGLKTCKISNSAEDDFGREPEYFKSRYGDISKLKRIYDASLNVEGLIV  
PDVSVHPWMQIKECSNYMYCAYATGDGHTNSELDEGRLITFTEKELRYINDISRTWVAN

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:71.2 kDa

Predicted isoelectric point:-

## 【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)


电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



**Storage:** Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

### 【Identification】



### 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

