

重组蛋白

Catalog Number:ORSBP4371R0

Recombinant Arabidopsis thaliana Protein NLP7(NLP7)

Organism Species: Arabidopsis thaliana (Mouse-ear
cress)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: NLP7

Uniprot:Q84TH9

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: Baculovirus

Protein Length: Full Length

Expression Region: 1-959aa

Sequence:

MCEPDDNSARNGVTTQPSRSRELLMDVDDLDDGSWPLDQIPYLSSSNRMISPIFVSSSS
EQPCSPLWAFSDGGNGFHHATSGGDDEKISSVSGVPSFRLAEYPLFLPYSSPSAAENTT
EKHNSFQFPSPLMSLVPPENTDNYCVIKERMTQALRYFKESTEQHVLAQVWAPVRKNGRD
LLTTLGQPFVLNPNGLNQRMISLTYMFSVDSESDVELGLPGRVFRQKLPEWTPNVQY
YSSKEFSRLDHALHYNVRGTLALPVFNPSGQSCIGVVELIMTSEKIHYAPEVDKVKCALE
AVNLKSSEILDHQTTQICNESRQNALAEILEVLTVCETHNLPLAQTWVPCQHGSVLANG
GGLKKNCTSFDGSCMGQICMSTTDMACYVVDAAHVWGFRDACLEHHLQKGQGVAGRAFLNG
GSCFCRDITKFCKTQYPLVHYALMFKLTTCAISLQSSYTGDDSYILEFFLPSSITDDQE
QDLLLGSILVTMKEHFQSLRVASGVDFGEDDDKLSFEIIQALPKKVVHASKIESIRVPFSG
FKSNATETMLIPQPVVQSSDPVNEKINVATVNGVVKEKKKTEKKRGKTEKTISLDVLQQY
FTGSLKDAKSLGVCPTTMKRICRQHGISRWPSRKIKKVNRSITKLKRVIESVQGTGGGL
DLTSMVSSIPWTHGQTSAPLNSPNGSKPELPNTNNSPNHWSSDHSPNEPNGSPELPP
SNGHKRSRTVDESAGTPTSHGSCDGNQLDEPKVPNQDPLFTVGGSPGLLFPYSRDHDVS
AASFAMPNRLLSIDHFRGMLIEDAGSSKDLRNLCPAAFDKQDTNWMNNDNNSNNNL
YAPPKEEAIANVACEPSGSEMRTVTIKASYKDDIIRFRISGSGIMELKDEVAKRLKVDA
GTFDIKYLDDDNEWVLIACDADLQECLEIPRSSRTKIVRLLVHDVTTNLGSSCESTGEL

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:107.4 kDa

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

