

重组蛋白

Catalog Number:ORSCF1421R0

Recombinant Pongo abelii NADPH oxidase 4(NOX4)

Organism Species: Pongo abelii (Sumatran orangutan)

(Pongo pygmaeus abelii)



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: NOX4

Uniprot:Q5R5C5

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-578aa

Sequence:

MAVSWRSWLANEGVKHLCLFIWLSMNVLFWKTFLLYNQGPYHYLHQMLGLGLCLSRAS
ASVLNLCNLSILLPMCRTLAYLRGSQKVPSRRTRRLDKSRTFHITCGVTICIFSGVHV
AAHLVNALNFSVNYSEDFVELNAARYRDEDPRKLLFTTVPGLTGVCMMVVFLMITASTY
AIRVSNYDIFWYTHNLFFVFYMLLTLHVSGLLKYQTNLDTHPPGCISLNRRTSSQNISLP
EYFSEHFHEPFPEGFSKPEEFTQNTFVKICMEEPFRQANFPQTLWISGPLCLYCAERLY
RYIRSNKPVTIISVISHPSDVMEIRMVKENFKARPGQYITLHCPSVSALENHPFTLT MCP
TETKATFGVHLKIVGDWTERFRDLLPPSSQDSEILPFIQSRNYPKLYIDGPF GSPFEES
LNYEVSLCVAGGIGVTPFASILNLLDDWKPYKLRRLYFIWVCRDIQSFRWFADLLCMLH
NKFWQENRPDYVNIQLYLSQTDGIQKIIGEKYHALNSRLFGRPRWKLLFDEIAKYNRGK
TVGVFCCGPNLSKTLHKLNSQINSYGTRFEYNKESFS

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:69.7 kDa

Predicted isoelectric point:-

【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

