

重组蛋白

Catalog Number:ORSCF2095RHu2

Recombinant Human Nucleolar transcription factor

1 (UBTF)

Organism Species: Homo sapiens (Human)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: UBTF

Uniprot:P17480

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-764aa

Sequence:

MNGEADCPTDLEMAAPKGQDRWSQEDMLTLLECMKNNLPSNDSSKFKTTESHMDWEKVAF
KDFSGDMCKLKWVEISNEVRKFRTLTELILDAQEHVKNPYKGGKLLKHPDFPKKPLTPYF
RFFMEKRAKYAKLHPEMSNDLTKILSKKYKELPEKKMKYIQDFQREKQEFERNLARFR
EDHPDLIQNAKSDIPEKPKTPQQLWYTHEKKVYLKVRPDATTKEVKDSLGGKQWSQLSDK
KRLKWIHKALEQRKEYEEMRDYIQKHPELNISEGITKSTLTKAERQLKDKFDGRPTKP
PPNSYSLYCAELMANMKDVPSTERMVLCSQQWKLSSQEKDAYHKKCDQKKKDYEVLLR
FLESLPEEEQQRVLGEEKMLNINKKQATSPASKKPAQEGGKGGSEKPKRPVSAMFIFSEE
KRRQLQEERPELSESELTRLLARMWNDLSEKKKAKYKAREAALKAQSERKPGGEREERGK
LPESPKRAEEIWQQSVIGDYLARFKNDRVKALKAMEMTWNMEKKEKLMWIKKAAEDQKR
YERELSEMRAPPAATNSSKKMKFQGEPPKPPMNGYQKFSQELLSNGELNHLPLKERMVEI
GSRWQRISQSQKEHYKLAEEQKQYKVHLDLWVKSLSPODRAAYKEYISNKRKSMTKLR
GPNPKSSRTTLQSKSESEEDDE
GDENEDE
DE

Tag Info: N-terminal 10xHis-SUMO-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:107.9 kDa

Predicted isoelectric point:-

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

