

重组蛋白

Catalog Number:ORSCF2457R0

Recombinant Mycoplasma pneumoniae P30 adhesin(p30)

Organism Species: Mycoplasma pneumoniae (strain ATCC
29342 / M129)



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: p30

Uniprot:P75330

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-274aa

Sequence:

MKLPPRRKCLKLFLLAWMLVLFSALIVLATLILVQHNTELTEVKSELSPLNVVLHAEEDT
VQIQGKPITEQAWFIPTVAGCFGFSALAIILGLAIGLPIVKRKEKRLLEEKERQEQLAEQ
LQRISAQQEEQQALEQQAAAEAHAEAEVEPAPQPVPVPPQPQVQINFGPRTGFPPQPGMA
PRPGMPPHPGMAPRPGFPPQPGMAPRPGMPPHPGMAPRPGFPPQPGMAPRPGMPPHPGMA
PRPGFPPQPGMAPRPGMQPPRPGMPPQPGFPPKR

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:32.6 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Identification】

kDa M



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

