

# 重组蛋白

Catalog Number:ORSCF3094R0

Recombinant Mycoplasma genitalium Uncharacterized  
protein MG281(MG281)

Organism Species: Mycoplasma genitalium (strain ATCC  
33530 / G-37 / NCTC 10195)



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

## 【Product Information】

Gene Names: MG281

Uniprot:P47523

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-556aa

Sequence:

MQFKKHKNSVKFKRKLFWTIGVLGAGALTFSAVMITNLVNQSGYALVASGRSGNLGFKL  
FSTQSPSAEVKLSLSLNDGSYQSEIDLSSGGANFREKFRNFANELSEAITNSPKGLDRPV  
PKTEISGLIKTGNFITPSFKAGYYDHVASDGSLLSYYQSTEYFNNRVLMPILQTTNGTL  
MANNRGYDDVFRQVPSFSGWSNTKATTVSTSNNLT YDKWTFYFAAKGSPLYDSYPNHFFED  
VKTLAIDAKDISALKTTIDSEKPTYLIIRGLSGNGSQLNELQLPESVKKVSLEYDYTGVN  
VAKQIFANVVELEFYSTSKANSFGFNPLVLGSKTNVIYDLFASKPFTHTIDLTQVTLQNSD  
NSAIDANKLKQAVGDIYNYRRFERQFQGYFAGGYIDKYLKKNVNTNKDSDDDLVYRSLKE  
LNLHLEEAYREGDNTYYRVNENYYPGASIEYENERASRDSEFQNEILKRAEQNGVTFDENI  
KRITASGKYSVQFQKLENDTDSLSLERTKAVEGLVTVIGEEKFETVDITGVSSDTNEVKS  
LAKELKTNALGVKLLK

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:65.0 kDa

Predicted isoelectric point:-

## 【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)


电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



**Storage:** Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

### 【Identification】



### 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

