

重组蛋白

Catalog Number:ORSCF4449RHu

Recombinant Human TIR domain-containing adapter
molecule 1(TICAM1)

Organism Species: Homo sapiens (Human)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: TICAM1

Uniprot:Q8IUC6

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-712aa

Sequence:

MACTGPSLPSAFDILGAAGQDKLLYLKHKLKTTPRPGCQGQDLLHAMVLLKLGQETEARIS
LEALKADAVARLVARQWAGVDSTEDPEEPPDVSWAVARLYHLLAEEKLCPASLRDVAIQE
AVRTLSSRDDHRLGELQDEARNRCGWDIAGDPGSIRTLQSNLGCPLPSSALPSGTRSLPR
PIDGVSDWSQGCSLRSTGSPASLASNLEISQSPTMPFLSLHRSPHGPSKLCDDPQASLVP
EPVPGGCQEPEEMSWPPSGEIASPPELPSSPPPGLPEVAPDATSTGLPDTAAPETSTNY
PVECTEGSAGPQSLPLPILEPVKNPCSVKDQTPQLQSVEDTTSPTNKPCPPTPTTPETSP
PPPPPPSSTPCSAHLTPSSLFPSSLESSEQKFYNFVILHARADEHIALRVREKLEALG
VPDGATFCEDFQVPGRGELSCLQDAIDHSAFIILLTNSNFCRLSLHQVNQAMMSNLTRQ
GSPDCVIPFLPLESSPAQLSSDTASLLSGLVRLDEHSQIFARKVANTFKPHRLQARKAMW
RKEQDTRALREQSQHLDGERMQAAALNAAAYSAYLQSYLSYQAQMEQLQVAFGSHMSFGTG
APYGARMPFGGQVPLGAPPPFPTWPGCPQPPPLHAWQAGTPPPSPQPAAFPQSLPFPQS
PAFPTASPAPPQSPGLQPLIHHAQMVQLGLNNHMWNQRGSQAPEDKTQEAE

Tag Info: N-terminal 10xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:79.9 kDa

Predicted isoelectric point:-

【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信