

重组蛋白

Catalog Number:ORSEP0289R0

Recombinant Dermatophagoides pteronyssinus Major mite

allergen Der p 23

Organism Species: Dermatophagoides pteronyssinus

(European house dust mite)

使用说明手册V1.0

 网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: N/A

Uniprot:L7N6F8

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length of Mature Protein

Expression Region: 22-90aa

Sequence:

ANDNDDDPTTTVHPTTTEQPDDKFECPSRFGYFADPKDPHKFYICSNWEAVHKDCPGNTR
WNEDEETCT

Tag Info: N-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:13.0 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Identification】

kDa M



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

