

重组蛋白

Catalog Number:ORSEP0299R0

Recombinant Arachis hypogaea Allergen Arachin Arah3

isoform

Organism Species: Arachis hypogaea (Peanut)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: N/A

Uniprot:B5TYU1

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length of Mature Protein

Expression Region: 21-530aa

Sequence:

ISFRQQPEENACQFQRLNAQRPDNRIESEGGYIETWNPNNQEFECAGVALSRLVLRNAL
RRPFYSNAPQEIFIQQGRGYFGLIFPGCPSTYEPAQQGRRYQSQRPPRRLQEEDQSQQQ
QDSHQKVHRFNEGDLIAVPTGVAFWLYNDHDTDVVAVSLTDTNNDNQLDQFPRRFNLAG
NHEQEFLRYQQQSRQSRRRSLPYSYSPQSQPRQEEREFSPRGQHSRRERAGQEEENEGG
NIFSGFTPEFLAQAFQVDDRQIVQNLRGENESEEQGAIVTVRGLRILSPDRKRGADEEE
EYDEDEYDEEDRRRGRGSRGSGNGIETICTATVKKNIGRNRSPDIYNPQAGSLKTAN
ELNLLILRWLGLSAEYGNLYRNALFVPHYNTNAHSIIYALRGRAHVQVVDSDNGNRVYDEE
LQEGHVLVVPQNFVAVAGKSQSDNFEYVAFKTDSPSIANLAGENSVIDNLPEEVVANSYG
LPREQARQLKNNNPFFKFFVPPSQQSPRAVA

Tag Info: N-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:62.5 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

