

重组蛋白

Catalog Number:ORSEP0410RHu

Recombinant Human CD44 antigen(CD44), partial

Organism Species: Homo sapiens (Human)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: CD44

Uniprot:P16070

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Partial of BC004372

Expression Region: 21-606aa

Sequence:

QIDLNITCRFAGVFHVEKNGRYSISRTEAADLCKAFNSTLPTMAQMEKALSIGFETCRYG
FIEGHVVIPRIHPNSICAAANTGVYILTSNTSQYDTYCFNASAPPEEDCTSVTDLPNAFD
GPITITIVNRDGTRYVQKGEYRTNPEDIYPSNPTDDDVSSEGSSERSSTSGGYIFYTFST
VHIPDEDSPWITDSTDRIPATSTSSNTISAGWEPNEENEDERDRHLSFSGSGIDDEDF
ISSTISTTPRAFDHTKQNQDWTQWNPSHSNPEVLLQTTTRMTDVDRNGTTAYEGNWNPEA
HPPLIHHEHHEEEETPHSTSTIQATPSSSTEETATQKEQWFGNRWHEGYRQTPREDSHST
TGTAASAHTSHPMQGRTPSPEDSSWTDFFENPISHPMGRGHQAGRRMDMDSSHSTTLQP
TANPNTGLVEDLDRTGPLSMTTQQSNSQSFSTSHGLEEDKDHPPTSTLTSSNRNDVTGG
RRDPNHSEGSTTLLEGYTSHPHTKESRTFIPVTSAKTGSFGVTAVTVGDSNSNVNRSLS
GDQDTFHPSGGSHHTHGSESDGHSHGSQEGGANTTSGPIRTPQIPE

Tag Info: N-terminal 6xHis-SUMO-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:80.2 kDa

Predicted isoelectric point:-

【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

