

# 重组蛋白

Catalog Number:ORSEP2466R0

Recombinant *Synechocystis* sp. Carbon dioxide-  
concentrating mechanism protein CcmK homolog 2(ccmK2)

Organism Species: *Synechocystis* sp. (strain PCC 6803 /  
Kazusa)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

**【Product Information】**

Gene Names: ccmK2

Uniprot:P72761

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length of Mature Protein

Expression Region: 2-103aa

Sequence:

SIAVGMIETRGFPVAVEAADSMVKAARVTLVGYEKIGSGRVTIVIRGDVSEVQASVSAGI  
EAANRVNGGEVLSTHIIARPHENLEYVLPPIRYTEEEVEQFRTY

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:18.4 kDa

Predicted isoelectric point:-

**【Storage and Stability】**

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)


邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



## 【Identification】



## 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

