

# 重组蛋白

Catalog Number:ORSEP3063R0

Recombinant Escherichia coli Chaperone protein

fimC(fimC)

Organism Species: Escherichia coli (strain K12)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

## 【Product Information】

Gene Names: fimC

Uniprot:P31697

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length of Mature Protein

Expression Region: 37-241aa

Sequence:

GVALGATRVIYPAGQKQEQ LAVTNNDENSTYLIQSWVENADGVKDGRFIVTPPLFAMK GK  
KENTLRILDATNNQLPQDRESLFWMNVKAIPSM DKS KLTENTLQLAIISRIKLYYRPAKL  
ALPPDQAAEKLRFRRSANSLTLINPTPYLLTVTELNAGTRVLENALVPPMGESTVKLPSD  
AGSNITYRTINDYGALTPKMTGVME

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:30.2 kDa

Predicted isoelectric point:-

## 【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)


邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



## 【Identification】



## 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

 网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

