

重组蛋白

Catalog Number:ORSEP3123AMu

Recombinant Escherichia coli Constitutive lysine

decarboxylase(ldcC) (Active)

Organism Species: Escherichia (strain K12)

使用说明手册V1.0

网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>

手机号: 13061929181(微信)
13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names:

Uniprot:P52095

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: Not test.

Source: E.coli

Protein Length: Full Length

Expression Region: 1-713aa

Sequence:

MNIIAIMGPHGVFYKDEPIKELESALVAQGFQIIWPQNSVDLLKFIEHNPRICGVIFDWD
EYSLDLCSDINQLNEYLPYAFINTHSTMDVSVQDMRMALWFFEYALGQAEDIAIRMRQY
TDEYLDNITPPFTKALFTYVKERKYTFCTPGHMGGTAYQKSPVGCLFYDFFGGNTLKADV
SISVTELGSLLDHTGPHLEAEEYIARTFGAEQSYIVTNGTSTSNKIVGMYAAPSGSTLLI
DRNCHKSLAHLMMNDVVPVWLKPTRNALGILGGIPRREFTRDSIEEKVAATTQAQWPVH
AVITNSTYDGLLYNTDWIKQTLDVPSIHFD SAWVPYTHFHP IYQKGKSGMSGERVAGKVIF
ETQSTHKMLAALSQASLIHIKGEYDEEAFNEAFMMHTTTSPSYPIVASVETAAAML RGNP
GKRLINRSVERALHFRKEVQRLREESDGWFFDIWQPPQVDEAECWPVAPGEQWHGFNDAD
ADHMF LDPVKVTILTPGMDEQGNMSEEGIPAALVAKFLDERGIVVEKTGPYNLLFLFSIG
IDKTKAMGLLRGLTEFKRSYDLNLRIKNMLPDLYAEDPDFYRNMRIQDLAQGIHKLIRKH
DLPGLMLRAFDTLPEMIMTPHQAWQRQIKGEVETIALEQLVGRV SANMILPYPPGVPLLM
PGEMLTKE SRTVLD FLLMLCSVGQHYPGFETDIHGAKQDEDGVYRVRVLKMAG

Tag Info: C-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.

If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Applications:-

Predicted Molecular Mass:-



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

Accurate Molecular Mass:82.6 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】

【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信