

重组蛋白

Catalog Number:ORSEP3570R0

Recombinant Mycobacterium bovis Diacylglycerol
acyltransferase/mycolyltransferase Ag85A(fbpa)

Organism Species: Mycobacterium bovis (strain BCG /
Pasteur 1173P2)

使用说明手册V1.0

 网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: fbpA

Uniprot:A1KQD8

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length of Mature Protein

Expression Region: 43-338aa

Sequence:

AFSRPGLPVEYLQVSPSPMGRDIKVQFQSGGANSPALYLLDGLRAQDDFSGWDINTPAFE
WYDQSGLSVVMPVGGQSSFYSDWYQPACGKAGCQTYKWETFLTSELPGWLQANRHVKPTG
SAVVGLSMAASSALTLAIYHPQQFVYAGAMSGLLDPSQAMGPTLIGLAMGDAGGYKASDM
WGPKEDPAWQRNDPLLNVGKLIANNTRVWVYCGNGKPSDLGGNNLPAKFLEGFVRTSNIK
FQDAYNAGGGHNGVDFDPDSGTHSWEYWGAQLNAMKPDQLRALGATPNTGPAPQGA

Tag Info: N-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:35.8 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Identification】

kDa M



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

