

重组蛋白

Catalog Number:ORSEP4076R0

Recombinant Staphylococcus epidermidis Shikimate

dehydrogenase (NADP(+)) (aroE)

Organism Species: Staphylococcus epidermidis (strain

ATCC 35984 / RP62A)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: aroE

Uniprot:Q5HNV1

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: E.coli

Protein Length: Full Length

Expression Region: 1-269aa

Sequence:

MKFAVIGNPISHSLSPMHANFQSLNLENTYEAINVPVNQFQDIKKIISEKSIDGFNVT
IPHKERIIPYLDDINEQAQSVGAVNTVLVKDGKWIGYNTDGIGYVNGLKQIYEGIEDAYI
LILGAGGASKGIANELYKIVRPTLTVANRTMSRFNNWSLNINKINLSHAESHLDEFDIII
NTTPAGMNGNTDSVISLNRLASHTLVSDIVYNPYKTPILIEAEQRCNPIYNGLDMFVHQQ
AESFKIWTNLEPDIKAMKNIVIQKLKGEL

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:37.4 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552


发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>

手机号: 13061929181(微信)
13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信