

重组蛋白

Catalog Number:ORSEP4455R0

Recombinant *Danio rerio* Interferon-induced GTP-binding

protein MxA (mxa), partial

Organism Species: *Danio rerio* (Zebrafish) (*Brachydanio rerio*)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: mxa

Uniprot:Q8JH68

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: partial

Expression Region: 31-310aa

Sequence:

VEKDLALPAIAVIGDQSSGKSSVLEALSGVPLPRGSGIVTRCPELEKLMIRTKDQDRWHGR
ISYKTCEEDFDDPAEVEKKIRQAQDEMAGAGVGISEELISLQITSADVPDLTLIDLPGIA
RVAVKGPENIGDQIKRLIRKFVTRQETINLVVPCNVDIATTEALQMAQAEDPDGERTL
GILTKPDLVDKGTEGTVVDIVHNEVIHLTKGYMIVRCRGQKEIMDQVTLNEATETESAFF
KDHPHFSKLYEEGFATIPKLAEKLTIELVHHIQKSLPRLE

Tag Info: N-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:35.0 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C. Avoid repeated

freeze/thaw cycles.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)


邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

