

重组蛋白

Catalog Number:ORSEP4522RHu3

Recombinant Human NACHT, LRR and PYD domains-containing
protein 3(NLRP3), partial

Organism Species: Homo sapiens (Human)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



【Product Information】

Gene Names: NLRP3

Uniprot:Q96P20

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: partial

Expression Region: 733-1036aa

Sequence:

LFSVLSTSQSLTELDLSDNSLGDPGMRVLCETLQHPGCNIRRLWLGRCGLSHECCFDISL
VLSSNQKLVELDLSNALGDFGIRLLCVGLKHLKLLCNLKKLWLVSCCLTSACQDLASVLS
TSHSLTRLYVGENALGDSGVAILCEKAKNPQC�NLQKLGLVNSGLTSVCCSALSSVLSTNQ
NLTHLYLRGNTLGDKGIKLLCEGLLHPDCKLQVLELDNCNLTSHCCWDLSTLLTSSQSLR
KLSLGNNDLGDVGMMFCEVLKQQSCLLQNLGLSEMYFNKYETKSALETLQEEKPELTVVF
EPSW

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:40.7 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

