

重组蛋白

Catalog Number:ORSEP4770R0

Recombinant Escherichia coli 0157:H7 ATP-dependent DNA

helicase recG(recG)

Organism Species: Escherichia coli 0157:H7



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: recG

Uniprot:Q8XD86

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length

Expression Region: 1-693aa

Sequence:

MTGRLLDAVPLSSLTGVGAALSNKLAKINLHTVQDLLLHLPLRYEDRTHLYPIGELLPGV
YATVEGEVLNLCNISFGGRRMMTCQISDGSILTMRFNFSAAMKNSLATGRRVLAYGEAK
RGKYGAEMIHPEYRVQGDLPSTPELQETLTPVYPTTEGVKQATLRKLTQALDLLDTCAIE
ELLPPELSQGMMLPEALRTLHRPPPTLQLSDLETGQHPAQRRLILEELLAHNLSMLALR
AGAQRFAQPLSANDALKNKLLAALPFKPTGAQARVVAEIERDMALDVPMMRLVQGDVGS
GKTLVAALAALRAIAHGKQVALMAPTELLAEQHANNFRNWFEPGLIEVGWLAGKQKGGKAR
LSQQEAIASGQVQMIVGTHAIFQEQVQFNGLALVIIDEQHRFGVHQRLALWEKGQQQGFH
PHQLIMTATPIPRTLAMTAYADLDTSVIDELPPGRTPVTTVAIPDTRRTDIIDRVRHACI
TEGRQAYWVCTLIEESELLEAQAAEATWEELKLALPELVGLVHGRMKPAEKQAVMASFK
QGELHLLVATTVIEVGVDVNPASLMI IENPERLGLAQLHQLRGRVGRGAVASHCVLLYKT
PLSKTAQIRLQVLRDSNDGFVIAQKDLEIRGPGELLGTRQTGNAEFKVADLLRDQAMIPE
VQRLARHIHERYPQQA KALIERWMPETERYSNA

Tag Info: N-terminal 6xHis-SUMO-tagged

Form: Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:89.4 kDa

Predicted isoelectric point:-

【Storage and Stability】

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



Storage: Store at 2-8°C for one month, the shelf life of liquid form is 6 months at -20° C/-80° C, The shelf life of lyophilized form is 12 months at -20° C/-80° C. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

