

# 重组蛋白

Catalog Number:ORSEP5008R0

Recombinant Arabidopsis thaliana UDP-  
glycosyltransferase 72B1 (UGT72B1)

Organism Species: Arabidopsis thaliana (Mouse-ear  
cress)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

 电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

## 【Product Information】

Gene Names: UGT72B1

Uniprot:Q9M156

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length

Expression Region: 1-480aa

Sequence:

MEESKTPHVAIIPSPGMGHLIPLVEFAKRLVHLHGLTVTFVIAGEGPPSKAQRTVLDSLPS  
SSISSVFLPPVDLTDLSSSTRIESRISLTVTRSNPELRKVFDSFVEGGRLPTALVVDLFG  
TDAFDVAVEFHVPPYIFYPTTANVLSFFLHLPKLDVTSCEFRELTEPLMLPGCVPVAGK  
DFLDPAQDRKDDAYKWLHNTKRYKEAEGILVNTFFELEPNAIKALQEPGLDKPPVYPVG  
PLVNIGKQEAQTEESECLKWLDNQPLGSVLYVSFGSGGTLTCEQLNELALGLADSEQRF  
LWVIRSPSGIANSSYFDSHSQTDPLTFLPPGFLERTKKRGFVIPFWAPQAQVLAHPSTGG  
FLTHCGWNSTLESVVSIGIPLIAWPLYAEQKMNAVLLSEDIRAALRPRAGDDGLVRREEVA  
RVVKGLMEGEEGKGVNRNKMKEAACRVLKDDGTSTKALSVALKWKAKHKELEQNGNH

Tag Info: N-terminal 10xHis-tagged and C-terminal Myc-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:60.4 kDa

Predicted isoelectric point:-

## 【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)


电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . Avoid repeated freeze/thaw cycles.

### 【Identification】



### 【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991

