

重组蛋白

Catalog Number:ORSEP5165R0

Recombinant Arabidopsis thaliana Pentatricopeptide

repeat-containing protein At3g42630 (At3g42630)

Organism Species: Arabidopsis thaliana (Mouse-ear
cress)



使用说明手册V1.0

网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

 邮箱: sale4@aoruisai.com
sale5@aoruisai.com

 电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: At3g42630

Uniprot:Q9M2A1

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: -

Source: E. coli

Protein Length: Full Length

Expression Region: 1-415aa

Sequence:

MLSLNLSLKPQHLKLLSCYTDSSAPSI AKKLIKESKLSRDFSQKIQIVDYAPLVQTL SQR
RLPDVAHEIFLQTKSVNLLPNYRTL CALMLCFAENGFVLRARTIWDEI INSCFVPDVFV
SKLISAYEQFGCFDEVAKITKDVAARHSKLLPVVSSLAISCFGKNGQLELMEGVIEEMDS
KGVLLAETANVIVRYYSFFGSLDKMEKAYGRVKKFGIVIEEEEIRAVVLAYLKQRKFYR
LREFLSDVGLGRRNLGNMLWNSVLLSYAADFKMKSLQREFIGMLDAGFSPDLTTFNIRAL
AFSRMALFWDLHLTLEHMRLNIVPDLVTFGCVVDAYMDKRLARNLEFVYNRMNLDDSP
VLTDPPLAFEVLGKGFHLSSEAVLEFSPRKNWTYRKLIGVYLKKKLRRDQIFWNY

Tag Info: N-terminal MBP-tagged and C-terminal 6xHis-tagged

Form: Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:91.5 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com


电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

